


About the two plots

Plot A: Original thresholds. Movement intensity never reaches the Walking threshold (0.020), so the system stays stuck in Queuing mode even while the walker is moving. Threshold was too high for our hardware.

Plot B: After lowering the threshold. System correctly switches from Queuing to Walking when the walker starts moving, and Beeping triggers when an obstacle is detected. Confirms the new threshold works.


Takeaway: Default thresholds need calibration per hardware — we tune them by observing actual sensor values during testing.

Same for distance threshold testing




Walker Proximity Sensor — Test Data

Distance to Nearest Obstacle (red vertical bars = buzzer fired)



Movement Intensity — above red = Walking, below orange = Queuing, yellow = no switch



System Decisions — which mode is active & when buzzer is armed

